




# MED15 Antibody

|                            |                                                                                                                                                                                                                                                                                                                                                                                                          |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA842753LA01HU                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                                                                                                                                                                                                                            |
| <b>Uniprot No.</b>         | Q96RN5                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Immunogen</b>           | Recombinant Human Mediator of RNA polymerase II transcription subunit 15 protein (433-530AA)                                                                                                                                                                                                                                                                                                             |
| <b>Raised In</b>           | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Species Reactivity</b>  | Human                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Tested Applications</b> | ELISA, WB, IF; Recommended dilution: WB:1:500-1:5000, IF:1:50-1:200                                                                                                                                                                                                                                                                                                                                      |
| <b>Form</b>                | Liquid                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Conjugate</b>           | Non-conjugated                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                                                                                                                                                                                                                                                                                                                         |
| <b>Purification Method</b> | >95%, Protein G purified                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Alias</b>               | Mediator of RNA polymerase II transcription subunit 15 (Activator-recruited cofactor 105 kDa component) (ARC105) (CTG repeat protein 7a) (Mediator complex subunit 15) (Positive cofactor 2 glutamine/Q-rich-associated protein) (PC2 glutamine/Q-rich-associated protein) (TPA-inducible gene 1 protein) (TIG-1) (Trinucleotide repeat-containing gene 7 protein), MED15, ARC105 CTG7A PCQAP TIG1 TNRC7 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Research Area</b>       | Epigenetics and Nuclear Signaling                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Target Names</b>        | MED15                                                                                                                                                                                                                                                                                                                                                                                                    |

## Image



### Western Blot

Positive WB detected in: MCF-7 whole cell lysate


All lanes: MED15 antibody at 5.2µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 87, 83, 76 kDa

Observed band size: 87 kDa



Immunofluorescence staining of SH-SY5Y cells with CSB-PA842753LA01HU at 1:100, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.