



TWNK Antibody, FITC conjugated

Product Code	CSB-PA842754LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RR1
Immunogen	Recombinant Human Twinkle protein, mitochondrial protein (559-684AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Twinkle protein, mitochondrial (EC 3.6.4.12) (Progressive external ophthalmoplegia 1 protein) (T7 gp4-like protein with intramitochondrial nucleoid localization) (T7-like mitochondrial DNA helicase) (Twinkle mtDNA helicase), TWNK, C10orf2 PEO1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TWNK