



NLGN4X Antibody, Biotin conjugated

Product Code	CSB-PA843131LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N0W4
Immunogen	Recombinant Human Neuroligin-4, X-linked protein (383-494AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuroligin-4, X-linked (Neuroligin X) (HNLX), NLGN4X, KIAA1260 NLGN4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NLGN4X
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.