



ZIC5 Antibody, Biotin conjugated

Product Code	CSB-PA850433LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96T25
Immunogen	Recombinant Human Zinc finger protein ZIC 5 protein (217-276AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein ZIC 5 (Zinc finger protein of the cerebellum 5), ZIC5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZIC5