



ZIC5 Antibody, Biotin conjugated

Product Code	CSB-PA850433LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96T25
Immunogen	Recombinant Human Zinc finger protein ZIC 5 protein (217-276AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein ZIC 5 (Zinc finger protein of the cerebellum 5), ZIC5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZIC5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.