



THEMIS Antibody, HRP conjugated

Product Code	CSB-PA850780LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N1K5
Immunogen	Recombinant Human Protein THEMIS protein (547-641AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein THEMIS (Thymocyte-expressed molecule involved in selection), THEMIS, C6orf190 C6orf207
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	THEMIS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.