



GBP4 Antibody, Biotin conjugated

Product Code	CSB-PA853481LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PP9
Immunogen	Recombinant Human Guanylate-binding protein 4 protein (494-603AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanylate-binding protein 4 (EC 3.6.5.-) (GTP-binding protein 4) (GBP-4) (Guanine nucleotide-binding protein 4), GBP4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GBP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.