



AQP10 Antibody, HRP conjugated

Product Code	CSB-PA853482LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PS8
Immunogen	Recombinant Human Aquaporin-10 protein (209-301AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Aquaporin-10 (AQP-10) (Aquaglyceroporin-10) (Small intestine aquaporin), AQP10
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	AQP10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.