



ALS2CR12 Antibody, HRP conjugated

Product Code	CSB-PA853485LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96Q35
Immunogen	Recombinant Human Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 12 protein (26-76AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 12 protein, ALS2CR12
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ALS2CR12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.