



EYA1 Antibody, HRP conjugated

Product Code	CSB-PA857862LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99502
Immunogen	Recombinant Human Eyes absent homolog 1 protein (168-321AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Eyes absent homolog 1 (EC 3.1.3.16) (EC 3.1.3.48), EYA1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	EYA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.