



SPAG8 Antibody, Biotin conjugated

Product Code	CSB-PA858730LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99932
Immunogen	Recombinant Human Sperm-associated antigen 8 protein (236-409AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sperm-associated antigen 8 (HSD-1) (Sperm membrane protein 1) (SMP-1) (Sperm membrane protein BS-84), SPAG8
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SPAG8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.