



IL12RB2 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA859518LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q99665 |
| Immunogen | Recombinant Human Interleukin-12 receptor subunit beta-2 protein (366-470AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Interleukin-12 receptor subunit beta-2 (IL-12 receptor subunit beta-2) (IL-12R subunit beta-2) (IL-12R-beta-2) (IL-12RB2), IL12RB2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | IL12RB2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |