



IL12RB2 Antibody, FITC conjugated

Product Code	CSB-PA859518LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99665
Immunogen	Recombinant Human Interleukin-12 receptor subunit beta-2 protein (366-470AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-12 receptor subunit beta-2 (IL-12 receptor subunit beta-2) (IL-12R subunit beta-2) (IL-12R-beta-2) (IL-12RB2), IL12RB2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL12RB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.