



COPS8 Antibody, Biotin conjugated

Product Code	CSB-PA860346LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99627
Immunogen	Recombinant Human COP9 signalosome complex subunit 8 protein (150-198AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	COP9 signalosome complex subunit 8 (SGN8) (Signalosome subunit 8) (COP9 homolog) (hCOP9) (JAB1-containing signalosome subunit 8), COPS8, CSN8
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	COPS8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.