



SLC12A5 Antibody, Biotin conjugated

Product Code	CSB-PA863965LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2X9
Immunogen	Recombinant Human Solute carrier family 12 member 5 protein (944-1071AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 12 member 5 (Electroneutral potassium-chloride cotransporter 2) (K-Cl cotransporter 2) (hKCC2) (Neuronal K-Cl cotransporter), SLC12A5, KCC2 KIAA1176
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC12A5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.