



SCUBE2 Antibody, Biotin conjugated

Product Code	CSB-PA865101LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQ36
Immunogen	Recombinant Human Signal peptide, CUB and EGF-like domain-containing protein 2 protein (534-703AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal peptide, CUB and EGF-like domain-containing protein 2 (Protein CEGP1) (Scube/You), SCUBE2, CEGP1
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	SCUBE2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.