



SEL1L Antibody, HRP conjugated

Product Code	CSB-PA866208LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBV2
Immunogen	Recombinant Human Protein sel-1 homolog 1 protein (22-181AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein sel-1 homolog 1 (Suppressor of lin-12-like protein 1) (Sel-1L), SEL1L, TSA305
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SEL1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.