



STRA6 Antibody, Biotin conjugated

Product Code	CSB-PA866301LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BX79
Immunogen	Recombinant Human Receptor for retinol uptake STRA6 protein (530-667AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor for retinol uptake STRA6 (Retinol-binding protein receptor STRA6) (Stimulated by retinoic acid gene 6 protein homolog), STRA6
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	STRA6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.