



SHARPIN Antibody, HRP conjugated

Product Code	CSB-PA867119LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0F6
Immunogen	Recombinant Human Sharpin protein (82-171AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sharpin (Shank-associated RH domain-interacting protein) (Shank-interacting protein-like 1) (hSIPL1), SHARPIN, SIPL1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SHARPIN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.