



ALG9 Antibody, Biotin conjugated

Product Code	CSB-PA867161LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H6U8
Immunogen	Recombinant Human Alpha-1,2-mannosyltransferase ALG9 protein (434-618AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-1,2-mannosyltransferase ALG9 (EC 2.4.1.259) (EC 2.4.1.261) (Asparagine-linked glycosylation protein 9 homolog) (Disrupted in bipolar disorder protein 1) (Dol-P-Man:Man(6)GlcNAc(2)-PP-Dol alpha-1,2-mannosyltransferase) (Dol-P-Man:Man(8)GlcNAc(2)-PP-Dol alpha-1,2-mannosyltransferase), ALG9, DIBD1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ALG9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.