



PNPLA8 Antibody, FITC conjugated

Product Code	CSB-PA868272LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NP80
Immunogen	Recombinant Human Calcium-independent phospholipase A2-gamma protein (158-345AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-independent phospholipase A2-gamma (EC 3.1.1.5) (Intracellular membrane-associated calcium-independent phospholipase A2 gamma) (iPLA2-gamma) (PNPLA-gamma) (Patatin-like phospholipase domain-containing protein 8) (iPLA2-2), PNPLA8, IPLA22 IPLA2G
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PNPLA8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.