



KIF15 Antibody, HRP conjugated

Product Code	CSB-PA868310LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NS87
Immunogen	Recombinant Human Kinesin-like protein KIF15 protein (1-129AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIF15 (Kinesin-like protein 2) (hKLP2) (Kinesin-like protein 7) (Serologically defined breast cancer antigen NY-BR-62), KIF15, KLP2 KNSL7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KIF15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.