



TRIM8 Antibody, Biotin conjugated

Product Code	CSB-PA871629LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZR9
Immunogen	Recombinant Human Probable E3 ubiquitin-protein ligase TRIM8 protein (367-517AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase TRIM8 (EC 2.3.2.27) (Glioblastoma-expressed RING finger protein) (RING finger protein 27) (RING-type E3 ubiquitin transferase TRIM8) (Tripartite motif-containing protein 8), TRIM8, GERP RNF27
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.