



TNKS1BP1 Antibody, Biotin conjugated

Product Code	CSB-PA871634LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C0C2
Immunogen	Recombinant Human 182 kDa tankyrase-1-binding protein (271-375AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	182 kDa tankyrase-1-binding protein, TNKS1BP1, KIAA1741 TAB182
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TNKS1BP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.