



MEGF9 Antibody, Biotin conjugated

Product Code	CSB-PA872434LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1U4
Immunogen	Recombinant Human Multiple epidermal growth factor-like domains protein 9 protein (49-178AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Multiple epidermal growth factor-like domains protein 9 (Multiple EGF-like domains protein 9) (Epidermal growth factor-like protein 5) (EGF-like protein 5), MEGF9, EGFL5 KIAA0818
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MEGF9