



KIF13B Antibody, Biotin conjugated

Product Code	CSB-PA873625LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQT8
Immunogen	Recombinant Human Kinesin-like protein KIF13B protein (1531-1701AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIF13B (Kinesin-like protein GAKIN), KIF13B, GAKIN KIAA0639
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KIF13B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.