



TNFRSF19 Antibody, HRP conjugated

Product Code	CSB-PA878878LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NS68
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 19 protein (276-415AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 19 (TRADE) (Toxicity and JNK inducer), TNFRSF19, TAJ TROY
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TNFRSF19
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.