



ONECUT1 Antibody, HRP conjugated

Product Code	CSB-PA880058LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBC0
Immunogen	Recombinant Human Hepatocyte nuclear factor 6 protein (192-292AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hepatocyte nuclear factor 6 (HNF-6) (One cut domain family member 1) (One cut homeobox 1), ONECUT1, HNF6 HNF6A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ONECUT1