



ONECUT1 Antibody, Biotin conjugated

Product Code	CSB-PA880058LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBC0
Immunogen	Recombinant Human Hepatocyte nuclear factor 6 protein (192-292AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hepatocyte nuclear factor 6 (HNF-6) (One cut domain family member 1) (One cut homeobox 1), ONECUT1, HNF6 HNF6A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ONECUT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.