



TXNIP Antibody, HRP conjugated

Product Code	CSB-PA880966LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H3M7
Immunogen	Recombinant Human Thioredoxin-interacting protein (299-381AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thioredoxin-interacting protein (Thioredoxin-binding protein 2) (Vitamin D3 up-regulated protein 1), TXNIP, VDUP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TXNIP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.