



EBF3 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA880977LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9H4W6 |
| Immunogen | Recombinant Human Transcription factor COE3 protein (416-530AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Transcription factor COE3 (Early B-cell factor 3) (EBF-3) (Olf-1/EBF-like 2) (O/E-2) (OE-2), EBF3, COE3 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | EBF3 |