



POPDC3 Antibody, Biotin conjugated

Product Code	CSB-PA881020LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBV1
Immunogen	Recombinant Human Popeye domain-containing protein 3 protein (130-291AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Popeye domain-containing protein 3 (Popeye protein 3), POPDC3, POP3
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	POPDC3