



CHMP1A Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA881032LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9HD42 |
| Immunogen | Recombinant Human Charged multivesicular body protein 1a protein (1-196AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Charged multivesicular body protein 1a (Chromatin-modifying protein 1a) (CHMP1a) (Vacuolar protein sorting-associated protein 46-1) (Vps46-1) (hVps46-1), CHMP1A, CHMP1 KIAA0047 PCOLN3 PRSM1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | CHMP1A |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |