



CHMP1A Antibody, FITC conjugated

Product Code	CSB-PA881032LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HD42
Immunogen	Recombinant Human Charged multivesicular body protein 1a protein (1-196AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Charged multivesicular body protein 1a (Chromatin-modifying protein 1a) (CHMP1a) (Vacuolar protein sorting-associated protein 46-1) (Vps46-1) (hVps46-1), CHMP1A, CHMP1 KIAA0047 PCOLN3 PRSM1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CHMP1A