



FANCL Antibody, Biotin conjugated

Product Code	CSB-PA882120LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NW38
Immunogen	Recombinant Human E3 ubiquitin-protein ligase FANCL protein (138-248AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase FANCL (EC 2.3.2.27) (Fanconi anemia group L protein) (Fanconi anemia-associated polypeptide of 43 kDa) (FAAP43) (RING-type E3 ubiquitin transferase FANCL), FANCL, PHF9
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FANCL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.