



BAZ2B Antibody, HRP conjugated

Product Code	CSB-PA883427LB01HU
Abbreviation	Bromodomain adjacent to zinc finger domain protein 2B
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIF8
Immunogen	Recombinant Human Bromodomain adjacent to zinc finger domain protein 2B protein (139-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	May play a role in transcriptional regulation interacting with ISWI.
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bromodomain adjacent to zinc finger domain protein 2B (hWALp4), BAZ2B, KIAA1476
Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	BAZ2B