



NANS Antibody, Biotin conjugated

Product Code	CSB-PA885699LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NR45
Immunogen	Recombinant Human Sialic acid synthase protein (132-259AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sialic acid synthase (N-acetylneuraminase synthase) (EC 2.5.1.56) (N-acetylneuraminase-9-phosphate synthase) (EC 2.5.1.57) (N-acetylneuraminic acid phosphate synthase) (N-acetylneuraminic acid synthase), NANS, SAS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NANS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.