



PNPO Antibody, Biotin conjugated

Product Code	CSB-PA885740LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NVS9
Immunogen	Recombinant Human Pyridoxine-5'-phosphate oxidase protein (114-207AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyridoxine-5'-phosphate oxidase (EC 1.4.3.5) (Pyridoxamine-phosphate oxidase), PNPO
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PNPO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.