




POLK Antibody

Product Code	CSB-PA887009LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBT6
Immunogen	Recombinant Human DNA polymerase kappa protein (227-337AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA polymerase kappa (EC 2.7.7.7) (DINB protein) (DINP), POLK, DINB1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POLK

Image



Western Blot
 Positive WB detected in: Hela whole cell lysate,
 Mouse kidney tissue
 All lanes: POLK antibody at 4.2µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 99, 55, 77, 43, 56, 53, 89
 kDa
 Observed band size: 53 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.