



USP24 Antibody, Biotin conjugated

Product Code	CSB-PA890777LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPU5
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 24 protein (2509-2620AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase 24 (EC 3.4.19.12) (Deubiquitinating enzyme 24) (Ubiquitin thioesterase 24) (Ubiquitin-specific-processing protease 24), USP24, KIAA1057
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	USP24
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.