



GJD2 Antibody, HRP conjugated

Product Code	CSB-PA891536LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKL4
Immunogen	Recombinant Human Gap junction delta-2 protein (99-197AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gap junction delta-2 protein (Connexin-36) (Cx36) (Gap junction alpha-9 protein), GJD2, GJA9
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GJD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.