



GJD2 Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA891536LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UKL4 |
| Immunogen | Recombinant Human Gap junction delta-2 protein (99-197AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Gap junction delta-2 protein (Connexin-36) (Cx36) (Gap junction alpha-9 protein), GJD2, GJA9 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | GJD2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |