



HAO1 Antibody, Biotin conjugated

Product Code	CSB-PA891938LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJM8
Immunogen	Recombinant Human Hydroxyacid oxidase 1 protein (161-273AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hydroxyacid oxidase 1 (HAOX1) (EC 1.1.3.15) (Glycolate oxidase) (GOX), HAO1, GOX1 HAOX1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	HAO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.