



TNRC6B Antibody, Biotin conjugated

Product Code	CSB-PA892489LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPQ9
Immunogen	Recombinant Human Trinucleotide repeat-containing gene 6B protein (639-779AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trinucleotide repeat-containing gene 6B protein, TNRC6B, KIAA1093
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TNRC6B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.