



# SHANK2 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA892492LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9UPX8
<b>Immunogen</b>	Recombinant Human SH3 and multiple ankyrin repeat domains protein 2 protein (728-822AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	SH3 and multiple ankyrin repeat domains protein 2 (Shank2) (Cortactin-binding protein 1) (CortBP1) (Proline-rich synapse-associated protein 1), SHANK2, CORTBP1 KIAA1022 PROSAP1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	SHANK2