



SHANK2 Antibody, FITC conjugated

Product Code	CSB-PA892492LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPX8
Immunogen	Recombinant Human SH3 and multiple ankyrin repeat domains protein 2 protein (728-822AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SH3 and multiple ankyrin repeat domains protein 2 (Shank2) (Cortactin-binding protein 1) (CortBP1) (Proline-rich synapse-associated protein 1), SHANK2, CORTBP1 KIAA1022 PROSAP1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SHANK2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.