



GTF3C5 Antibody, FITC conjugated

Product Code	CSB-PA896534LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5Q8
Immunogen	Recombinant Human General transcription factor 3C polypeptide 5 protein (305-448AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	General transcription factor 3C polypeptide 5 (TF3C-epsilon) (Transcription factor III C 63 kDa subunit) (TFIIIC 63 kDa subunit) (TFIIIC63) (Transcription factor III C subunit epsilon), GTF3C5
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	GTF3C5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.