



WASF2 Antibody, HRP conjugated

Product Code	CSB-PA896553LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6W5
Immunogen	Recombinant Human Wiskott-Aldrich syndrome protein family member 2 protein (243-317AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Wiskott-Aldrich syndrome protein family member 2 (WASP family protein member 2) (Protein WAVE-2) (Verprolin homology domain-containing protein 2), WASF2, WAVE2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WASF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.