



RAD54L2 Antibody, HRP conjugated

Product Code	CSB-PA896892LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4B4
Immunogen	Recombinant Human Helicase ARIP4 protein (1122-1296AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Helicase ARIP4 (EC 3.6.4.12) (Androgen receptor-interacting protein 4) (RAD54-like protein 2), RAD54L2, ARIP4 KIAA0809
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RAD54L2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.