



LAS1L Antibody, Biotin conjugated

Product Code	CSB-PA897285LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4W2
Immunogen	Recombinant Human Ribosomal biogenesis protein LAS1L protein (490-565AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribosomal biogenesis protein LAS1L (Protein LAS1 homolog), LAS1L
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LAS1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.